

Dry Ice solution just-in-time

Pharmaceutical products and biomedical specimens – especially body fluids taken for analysis during clinical trials – often require the use of dry ice as a coolant.

How it works

Investigators participating in a clinical trial order the dry ice they need by contacting the TNT Healthcare Centre I in their country. The dry ice arrives the next day before noon. The sample is then packed with the dry ice according to instructions, and TNT collects the shipment for delivery to the central laboratory in the afternoon of the same day. Our track & trace system facilitates monitoring of the shipment throughout its transit.

Using our Just-In-Time Dry Ice Solution eliminates the need to make arrangements with multiple suppliers with differing minimum order quantities, service levels and invoicing procedures.



Benefits

- Late cut-off time
- One centrally-managed service for most European countries
- A clear invoicing process
- Improved administrative control
- Reduced management time

TNT Healthcare Centre staff are trained to support your organisation in the distribution of biomedical samples, drugs and medical equipment.

Your shipment arrives in perfect condition.

